



# CATALOGUE



## WEEKLY PILL DISPENSING SYSTEM

- BLISTERS & LABELS
- ACCESSORIES
- DISTRIBUTION
- UNIDOSE



## A DAILY **PEACE OF MIND**

4 TAKEN PER DAY  
7-DAY DOSES

### ALL-IN-ONE

THE MOST EFFECTIVE  
MULTIDOSE DETACHABLE  
SYSTEM\*



- **Presentation:** Complete patient pharmacological information
- **Anchoring system:** Use positioning pegs to align the blister pack - simply the best
- Highly popular 8.5" x 11" format
- **Detachable**
- **Visual inspection:** transparent blister to view medications at a glance
- **Integrity:** medication protection assured by our effective labels
- **Prescription changes:** the pill dispenser can be modified and securely resealed
- Filled by hand or with automated system
- Braille-coded blisters
- **DISPILL Solutions:** offered in different sizes
- Treatment: 7 days / 4 daily doses per card
- Labels are colour-coded for easy recognition and use

The picture represents a fictional patient file

\*Compared with pouches and other blister packs



Eco-responsible  
Recyclable & shreddable



Made in Canada

# DISPILL HYBRID PERMANENT

DISPILL permanent **HYBRID** multidose system offers the patient

**TWO** simple ways to open.

**BEND**  
Row, Up & Down



**AND**

**DETACH**  
Row, then blister



**PUNCTURE**  
Along dotted line  
with your finger



**OR**

**PEEL**  
the label




**PULL the label to OPEN**



A full blister is easier to open  
using the movements **PULL** and **PEEL**!


*As easy to peel  
as a **tangerine**!*

A safe solution for patients with tremor or synchronicity problems with their hands.



## UNIQUE FEATURES

### Labels



#### Patented System

3-peg mounting method for perfect alignment of the label and blister card.



#### Anti-burr Edge

Produced with an anti-burr edge that controls adhesive residue when printing.

**High Performance Adhesive,**  
Label with a high performance permanent adhesive that adheres better to the blister.



#### Premium Fiber Labels

Paper produced from premium fiber for superior printing and greater stability, at any time.


**Brighter colours that enhance the time of dose and facilitates reading.**

### Blister cards

The right capacity for all your needs



Unidose Permanent Resealable Permanent Deeper Hypso Hybrid



**LARGEST CAPACITY  
ON THE MARKET**



**Blister cards with a precut edge**  
for ease of opening

### Automation




**DISPILL** unique features make it by far the best choice for pill dispensing automation.




## Hybrid Solution


300



4-11



**HYBRID LABEL**  
Permanent adhesive (300)  
or resealable adhesive (4-11)  
500 UN. / BOX



How to peel off the liner

300 Permanent Hybrid Label 11" L


4-11 Resealable Hybrid Label 11" L

## Permanent Solution

31-10




31-11



**PERMANENT LABEL**  
Tamper-evident  
Permanent seal

500 UN. / BOX



How to peel off the liner

31-10 Permanent Label 10" L

31-11 Permanent Label 11" L

250

250-P



**HYBRID BLISTER CARD**  
Clear  
For permanent or resealable  
hybrid label  
500 UN. / BOX

250 Hybrid Blister Regular

250-P Hybrid Blister Deeper

1

1-P



**PERMANENT BLISTER CARD**  
Clear  
For label with permanent adhesive  
500 UN. / BOX

1 Permanent Blister Regular



1-P Permanent Blister Deeper




**DISPILL PILL DISPENSING SYSTEM  
IS ECO-RESPONSIBLE**

recyclable • shreddable

## Resealable Solution



**RESEALABLE LABEL**  
With resealable adhesive  
for changes and cold  
pressure seal  
**500 UN. / BOX**





How to peel off the liner


**4-10 Resealable Label 10'' L**

**4-11 Resealable Label 11'' L**

## Hypso Solution



**METRIC FORMAT LABEL**  
210 X 297 mm  
Tamper-evident Permanent  
Adhesive with Hybrid  
opening (350) or with  
Resealable Adhesive (450)  
**500 UN. / BOX**



How to peel off the liner

**350 Permanent Hybrid Label 297 mm H**

**450 Resealable Label 297 mm H**



**RESEALABLE BLISTER CARD**  
Clear  
For label with resealable adhesive  
**500 UN. / BOX**

**201 Resealable Blister Regular**

**201-P Resealable Blister Deeper**




**HYBRID BLISTER CARD**  
Clear  
For permanent or resealable Hybrid label  
– Metric Format 210 x 225 mm  
**500 UN. / BOX**

**280 Hypso Hybrid Blister Deeper**

Note: Hypso Solution with its metric format requires specific mounting, verification and distribution accessories. Ask for them at your DISPILL Customer Service.

## Specialty Solution

41-10



**NARCOTIC LABEL**  
Permanent adhesive  
Tamper-evident  
Permanent seal  
100 OR 500 UN. / BOX

**41-10 Permanent Label for Narcotics 10" L**

**41-11 Permanent Label for Narcotics 11" L**

**42-10 Permanent Label for PRN 10" L**

**42-11 Permanent Label for PRN 11" L**

42-11




**PRN LABEL**  
Permanent adhesive  
Take as needed  
100 OR 500 UN. / BOX




## Unidose Solution

620



**CARDBOARD – 30 DAYS**  
Permanent adhesive  
With Foil membrane  
6" x 9" folded  
250 UN. / BOX




Sealed Unidose

**620 Cardboard Permanent – 30 days**

201-P-A

1-P-A

250-P-A



**AMBER BLISTER CARD**  
Translucent  
Amber blisters provide 100%  
UV protection to medications  
100 OR 500 UN. / BOX

**201-P-A Resealable Deeper Amber Blister**

**1-P-A Permanent Deeper Amber Blister**

**250-P-A Hybrid Deeper Amber Blister**

650



**BLISTER CARD**  
Clear  
For Cardboard  
Cold Seal  
250 UN. / BOX

**650 Unidose Blister Large – 30 days**

## PREPARATION

### SOLUTION COLD SEAL



The DISPILL multidose or unidose system is designed to instantly seal by compression. DISPILL requires no heat!

Creating a lightweight, durable and reliable solution. Our labels are made with specific adhesives to meet pharmaceutical requirements and ensure the integrity of the container throughout its use.

700



750 Hypso Mounting Tray

### 700 UNIVERSAL MOUNTING TRAY

Convenient for all DISPILL blister cards  
3 peg anchoring system  
White plastic

24



### MOUNTING TRAY SUPPORT STAND

Holds mounting tray at a perfect 3" angle  
For Universal or Hypso mounting tray  
Nickel finish

15



Apply constant pressure  
on the DISPILL by vertical  
movements,  
up and down.

### SEALING ROLLER

For cold pressure seal  
All labels  
Comfortable plastic handle  
Hard rubber roller

680



### UNIDOSE MOUNTING TRAY

Triple Cork mounting tray for 30/31/32 days  
Cork to be used for cold seal only



## VERIFICATION

19



### VERIFICATION TRAY

With sliding wings  
In white and black plastic

14



### VERIFICATION STAND

Can stack 13 to 15 blister cards  
For use before cold sealing  
Nickel finish

## CHANGES

50



### SELF-ADHESIVE SECURITY SEAL

Transparent sticker with a red alteration line  
To seal blister changes  
816 SEALS / PACK



37



### PRESSURE HAND TOOL

To strengthen the seal of the resealable label  
after a change  
Wood 6" soft and flat end  
3 UN. / PACK

35



### PRECISION SCALPEL

For precision cuts on a soft blister  
Stainless steel blade, plastic handle  
With retractable safety cap  
5 UN./ PACK


**Note :** To make a change safely, it is important to use the tools suggested by DISPILL and follow the recommended guidelines.



## From the pharmacy to patients



17



### STORAGE SUPPORT

For multifunction storage box  
In steel  
Can carry 28 to 36 complete DISPILL

16



### MULTIFUNCTION STORAGE BOX

Made of highly durable blue plastic. Shatter resistant interlocking lid. Option for security ties. Stackable. Format 49L 21" x 15" x 12"

39



### SECURITY TIES

To secure multifunction box  
In resistant nylon  
5,5"  
100 UN. / PACK

## COMMUNICATION

30



### DISPILL PRESCRIPTION BAG

Made of natural recyclable paper  
11" x 14"

500 UN. / PACK



**CAN HOLD  
4 DISPILL  
STACKED FLAT**



908



Instruction sheet  
pad - DISPILL Hybrid  
Full colour  
4" x 5"  
100 sheets per pad

## In nursing homes



13

### HOOP ROD STAND

For counter or wallmount  
Nickel finish  
Can hold 12 to 14 blisters  
10,5" H



System to flip through the DISPILL cards and tear off the required row or time of dose without ever having to remove the card from the stand. Each DISPILL hangs from a Polysheet Lifter.

108



### SEPARATORS

To be used with the hoop rod stand  
Plastified cardboard  
6 distinctive colours:  
Blue, white, lilac, pink, orange, green  
Sold separately


25



### POLY SHEET LIFTER

Needed to attach DISPILL cards to hoop rod.  
Clear plastic  
Sold separately

18



Attaches to the suspension bars of the drawer

Fits most standard size lateral filing drawers


### STAND FOR FILING CABINET

For DISPILL blister card  
Nickel finish  
3 rows of 15 to 18 DISPILL  
(45 to 54 DISPILL per drawer)

12-S

### TABLE TOP FILING STAND

For hanging DISPILL  
Nickel finish  
10,5" H x 8,5" L



12-S With dividers - holds 12 DISPILL

12 Holds 13-15 DISPILL

6



### DISTRIBUTION TRAY

In rigid white plastic  
35 positions for blisters  
7 positions for other medications

930

### MEDICATION CART

Dual lock model  
For storage and distribution of the DISPILL cards and all other medications drawers.

Removable dividers  
4 swivel castors

*Additional information and specifications on request  
Photo for reference only*



## OUR COMMITMENT



### PAYMENT

You can pay by cheque, PAP (pre-authorized payment) or EFT (Electronic Funds Transfer).

Terms of payment Net 30 days.

All of our products are taxable according to provincial laws.



### TECHNICAL EXPERTISE

We are available Monday to Friday.

Do not hesitate to contact us for any request, we will get back to you promptly.

Our team of representatives is always available to meet your needs.



### OUR SERVICES

- After sale service
- Product demos
- Staff training
- Order follow-up
- Inventory monitoring assistance
- Solution expert



### DELIVERY

Delivery options to suit your needs. From Monday to Friday.

Orders received before 2 p.m. ET will be shipped same day.

Express delivery where service is offered. Delivery charges may apply, ask your customer service rep.



### SATISFACTION GUARANTEED

Our Quality Control Managers rely on official standards and certifications to provide you with DISPILL products of the highest quality.

We will exchange or give you credit on defective items or shipping error. Certain conditions apply.



### ECO-RESPONSIBLE

All DISPILL products are manufactured taking into account environmental standards. DISPILL blisters and their labels are recyclable and shreddable. Transport boxes are also recyclable.

GO  
FOR IT!

**GODISPILL.COM**

## Online training made easy!

All our preparation methods and how to use DISPILL accessories are presented in short videos.



DISPILL, Method Modification Hybrid label | 1:30 min



Train faster  
and stay a  
step ahead





TO ORDER OR TO CONTACT US

Monday through Friday from 8:30 a.m. to 4:30 p.m. ET.



**Toll Free**  
**Granby**

1 877 347.7455  
450 372.1299



[infos@dispill.com](mailto:infos@dispill.com)



24 hours / 7 days



[dispill.com](http://dispill.com)

 **DISPILL®**

**DISPILL**

1008 Marcoux Street | Unit 154 | Granby, QC J2J 2M8

AN INNOVATION FROM



PHARMA SOLUTIONS